ri-magic® HPLED

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Instructions

Mode d'emploi

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Table of Contents



1. Introduction
- 1.1. Important information to note prior to initial setup.
- 1.2. Safety symbols
- 1.3. Packaging symbols
- 1.4. Intended use
- 1.4.1 Indication
- 1.4.2 Contraindication
- 1.4.3 Intended patient population
- 1.4.4 Intended operators/users
- 1.4.5 Required skills/operator training
- 1.4.6 Environmental conditions
- 1.5. Warnings/caution
- 1.6. User responsibility
- 1.7. Scope of delivery
2. Initial setup of the device
- 2.1. Device function
- 2.2. Purpose
3. Assembly instructions before initial setup
- 3.1. ri-magic®HPLED table model
- 3.2. ri-magic®HPLED wall model
- 3.3. ri-magic® HPLED floor model
- 3.4. Initial setup and function
4. Note on the flexible arm
5. Care instructions
6. Cleaning and disinfection
- 6.1. Accessories
7. Technical specifications
8. Maintenance
9. Standards
10. Disposal of packaging
11. Explosion hazard
12. Disposal of accessories and device
13. Accompanying documents on electromagnetic compatibility according to IEC 60601-1-2, 2014, Ed. 4.0
- 13.1. Disposal
14. Warranty

1. Introduction

1.1 Important Information Prior to Initial Setup.

You have purchased a high-quality Riester product, which has been manufactured in compliance with regulation (EU) 2017/745 and is always subject to the strictest quality controls. Read these instructions for use (IFU) carefully before using the device and keep them in a safe place. If you have any questions, we are available at any time, and our contact information is provided at the end of this IFU. The address of our sales partner can be obtained up on request. Please note all instruments described in these instructions for use may only be used by appropriately trained personnel. The safe functioning of this device is only guaranteed if Riester original parts and accessories are used.

1.2. Safety symbols

Symbol	Note on symbol
	Follow the instructions in the user manual.
	Type B applied part
	Medical device
	Class II protective devices
	Warning! The general warning symbol indicates a potentially dangerous situation that can lead to serious injuries.
	Caution! Important note in this manual. The caution symbol indicates a potentially dangerous situation that can lead to minor or moderate injuries. It can also be used to warn of unsafe practices.
	Caution: Do not look into the beam
	LED light Do not look into the beam Class 2 LED
	Do not use outdoors
	Direct current
	Alternating current
	Manufacturing date YYMMDD (year, month, day)
	Manufacturer
	Manufacturer's serial number
	Lot / batch number
	Reference number
	Temperature for transportation and storage
	Relative humidity for transportation and storage
	Air pressure for transportation and storage Operating ambient air pressure
	CE-marking
	Symbol for the marking of electrical and electronic equipment in accordance with Directive 2002/96/EC. Caution: Used electrical and electronic equipment should not be treated as normal household waste, but should be disposed of separately in accordance with national and EU regulations.
	Non-ionizing radiation
	Caution: (US) federal law restricts this device from being used by or by order of a doctor (licensed physician).

1.3 Packaging symbols

Symbol	Note on symbol
	Indicates that the contents of the transport package are fragile and should therefore be handled with care.
	Store in a dry place
	Indicates the correct position for transporting the package.
	Keep away from sunlight
	"Green Dot" (country-specific)

Notes on electromagnetic compatibility

There are currently no indications that electromagnetic interactions with other devices can occur when the devices are used as intended.

Nevertheless, under the increased influence of unfavourable field strengths, e.g. during the operation of wireless telephones or radiological instruments, disturbances cannot be completely excluded.

Warning:

The use of other accessories may result in increased electromagnetic emissions or reduced electromagnetic immunity of the device and lead to incorrect operation.

The electromagnetic compatibility of this device has been verified by tests according to the requirements of IEC 60601-1-2.

Caution!

Do not use the examination light for eye exams.

There may be a risk of ignition if the device is operated in the presence of flammable mixtures of pharmaceuticals and air, oxygen, nitrous oxide or anaesthetic gases!

The device may only be opened by persons with the appropriate specialist knowledge. There is a danger of life-threatening electrical shock.

Only use the plug-in power supply adapter supplied.

Manufacturer Glob Tek Type GTM96180-1507-2.0

1.4 Intended Use

The universally applicable ri-magic® HPLED examination lamp from Riester was made for lighting and thus for diagnostics in all areas of medicine. It provides illumination of >200,000 lux at a distance of 200 mm.

1.4.1 Indication

The ri-magic® HPLED examination lamp serves as a light source for diagnostic purposes in medical practices.

The various versions in the catalogue differ in mounting options. However, the functional electrical and housing structure are identical for all models.

Adaptation to different mains supply networks available, and easily engaged by the wide-range transformer.

The instrument aids doctors or specialists in the detection, diagnosis, monitoring, treatment or alleviation of illnesses, injuries, and disabilities.

1.4.2 Contraindication

There may be a risk of ignition of gases if the instrument is operated in the presence of flammable mixtures or mixtures of pharmaceuticals.

The examination lamps must never be placed in liquids.

Only use Riester or Riester-approved accessories/consumables.

Cleaning frequency and sequence must comply with the cleaning regulations of non-sterile products in the respective facility. Cleaning/disinfection instructions must be observed.

The product may only be used by trained personnel.

1.4.3 Intended patient population

The device is intended for use on adults and children.

1.4.4 Intended operators/users

The examination lamps are intended exclusively for use by doctors in clinics and medical practices.

1.4.5 Required skills/operator training

Only doctors / specialists / trained personnel should use the examination lamps, as they have the appropriate qualifications.

1.4.6 Environmental conditions

The device is intended for use in rooms with a controlled environment.
The device must not be exposed to adverse/harsh environmental conditions.

1.5. Warnings/caution

Warning!

The general warning symbol indicates a potentially dangerous situation that can lead to serious injuries.

Do not use in a magnetic resonance environment.

There is a risk of ignition of gases if the device is operated in the presence of flammable mixtures or mixtures of medicinal products and air or oxygen or nitrous oxide!
The device must not be operated in rooms in which flammable mixtures or mixtures of pharmaceuticals and air or oxygen or nitrous oxide are present, e.g. operating rooms.

Electric shock!
The housing of the examination lights may only be opened by authorized persons.

Damage to the device due to a fall or strong ESD influence!
If the device is not functioning, it must be returned to the manufacturer for repair.

The device must be used in a controlled environment.
The device must not be exposed to harsh environmental conditions.

Caution!

The caution symbol indicates a potentially dangerous situation that can lead to minor or moderate injuries. It can also be used to warn of unsafe practices.

The faultless and safe functioning of these examination lamps can only be guaranteed if Riester original parts and accessories are used.

Old electronic devices must be disposed of in accordance with the institutional guidelines for the disposal of expired devices.

Cleaning frequency and sequence must comply with the cleaning regulations of non-sterile products in the respective facility. Cleaning/disinfection instructions in the instructions for use must be observed.

We recommend unplugging the power supply before cleaning or disinfecting.
Clean and disinfect the examination lamps carefully so that no liquid penetrates the interior.
Never place the examination lamps in liquids!
The examination lamps are delivered in a non-sterile condition. Do not use ethylene, oxide gas, heat, autoclaves or any other harsh methods to sterilize the device. The devices have not been approved for mechanical reprocessing or sterilization. This leads to irreparable damage!

The patient is not the intended operator.
The product may only be used by qualified personnel.
Qualified personnel are doctors or nurses in hospitals, medical facilities, clinics and medical practices.

1.6. User Responsibility

Caution!

User Responsibility
It is your responsibility to:
The user must check the integrity and completeness of the examination lamps before each use. All components must be compatible with each other.

Incompatible components can result in degraded performance.
Never knowingly use a defective device.
Replace parts that are defective, worn, missing or incomplete.
Contact the nearest factory-approved service centre if repairs or replacements are required.
In addition, the user of the device bears sole responsibility for malfunctions resulting from improper use, incorrect maintenance, improper repair, damage or changes by persons other than Riester employees or authorized service personnel.

All serious incidents related to the product must be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

If there are any problems with the product or use of the product, please contact your supervisor immediately.

1.7. Scope of delivery

Article no.: 6230
Ri-magic HPLED, table model, including universal clamps 230 V/120

Article no.: 6230.001
Ri-magic HPLED, table model, including universal clamps 230 V/120 V with 3-pin plug

Article no.: 6230.002
Ri-magic HPLED, table model, including universal clamps 230 V/120 V with US plug

Article no.: 6230.003
Ri-magic HPLED, table model, including universal clamps 230 V/120 V with Australian plug

Article no.: 6240
Ri-magic HPLED wall model 230 V/120 V

Article no.: 6240.001
Ri-magic HPLED wall model 230 V/120 V with 3-pin plug

Article no.: 6240.002
Ri-magic HPLED wall model 230 V/120 V with US plug

Article no.: 6240.003
Ri-magic HPLED wall model 230 V/120 V with Australian plug

Article no.: 6250
Ri-magic HPLED floor model 230 V/120 V

Article no.: 6250.001
Ri-magic HPLED floor model 230 V/120 V with 3-pin plug

Article no.: 6250.002
Ri-magic HPLED floor model 230 V/120 V with US plug

Article no.: 6250.003
Ri-magic HPLED floor model 230 V/120 V with Australian plug

Article no.: 6252
Ri-magic HPLED floor model 230 V/120 V, power supply attached to stand

Article no.: 6252.001
Ri-magic HPLED floor model 230 V/120 V, British, power supply attached to stand

Article no.: 6252.002
Ri-magic HPLED floor model 230 V/120 V, Australian, power supply attached to stand

Article no.: 6252.003
Ri-magic HPLED floor model 230 V/120 V, US, power supply attached to stand

2. Initial setup of the device

2.1. Device function:

1. LED with magnifying lens
2. Focusing wheel
3. Cooling fins/body
4. On-off switch/rheostat for stepless adjustment of the light intensity


1. Flexible lamp arm
2. Screw clamp for table mounting
3. Locking screw

1. Wall mount
2. Jack plug - power supply connection

1. Cable with jack plug from the power supply
2. Bracket with power supply
3. Stand
4. Mobile stand
5. Castors

Lamp

Table model

Wall model

Ri-magic HPLED 6252 ri-magic HPLED 6250

2.2. Purpose

The universally applicable ri-magic® HPLED examination lamp from Riester was made for lighting and thus for diagnostics in all areas of medicine. It provides illumination of approximately 130,000 lux at a distance of 400 mm.

3. Assembly instructions before initial setup

3.1. ri-magic®HPLED table model

Attach the universal clamp (1) to the table or wall rail using the adjusting screw (2). Insert the pin contact of the lamp (4) into the hole in the universal clamp from above and tighten it with the knurled screw (5). Insert the cable connector (6) of the power supply (7) into the socket (8) from below until you feel it click into place. Plug the transformer (7) into the socket. The lamp is now functional.

3.2. ri-magic®HPLED wall model

Hold the wall mount (1a) in the desired location on the wall and mark 4 holes for the dowels / anchors. Remove the wall mount and drill the 4 marked holes. Insert dowels / anchors into the wall. Attach the wall mount and tighten. Insert the lamp (4) into the mount from above and it in place with the locking screw (3a) on the side. Insert the cable connector (6a) of the power supply (7a) into the socket (8a) from below until you feel it click into place. Plug the transformer (7a) into the socket. The lamp is now functional.

Caution:

The wall must be sufficiently stable to ensure that the examination light is securely attached to the wall. Please only use the screws and dowels / anchors supplied. Alternatively, 6 mm dowels / anchors can be used in conjunction with rounded head screws 3 x 40 mm DIN 95. When possible it is best to mount the exam light into a wall stud (as opposed to a hollow wall) for maximum durability.

Installation instructions:

Please pay attention to the enclosed assembly instructions of the mounting systems when installing the accessories.

In order to ensure complete disconnection from the power supply at all times, the ME device must be installed in such a way that the power cable is accessible and can be disconnected.

3.3 ri-magic® HPLED floor model

Follow the assembly instructions enclosed with the mobile base to assemble the mobile base. Insert the cable connector (6b) of the power supply (7b) into the socket (8b) from the side until you hear a click. The lamp is now functional.

To be able to move the light from one location to another, the flexible arm (head) must be bent downwards, so that the light is not too high and does not hit nearby objects and so that there is no excess weight when the light is slightly inclined.

3.4. Initial setup and function

First insert the cable connector (6a/6b) into the socket (8a/8b). Connect the power supply (7a/7b) to the mains via the socket.

On/off

Turn the rheostat (9) clockwise. When you hear a click, the light is switched on. Turn the rheostat (9) – counter-clockwise until you hear a click. The light is switched off. The examination light can also be disconnected from the mains by removing the power adapter from the socket.

Rheostat

You can set the desired light intensity by turning the rheostat (9). Turning it clockwise increases the light intensity, turning it counter-clockwise decreases it.

Focusing

A focusing ring (10) is located on the front of the lamp head. This ring can be turned to adjust the coverage area of the light.

4. Note on the flexible arm

The flexible lamp arm should not be bent too tightly, as this could lead to premature material fatigue.

During normal operation of the light, the mount (4) and the universal clamp (1) or stand (17) will become warm in the area of the socket (8a/8b).

5. Care instructions

General note

The cleaning and disinfecting of the medical devices serve to protect the patient, the user and third parties and to maintain the value of the medical devices. Due to the product design and the materials used, a defined limit for the maximum possible number of re-processing cycles cannot be determined. The service life of medical devices is defined by their function and careful handling. Before return for repair, defective products must have undergone the prescribed reconditioning process.

6. Cleaning and disinfection

In order to avoid possible cross-contamination, the examination lamps must be cleaned and disinfected regularly.

The examination lamps can be cleaned on the outside using a damp cloth (if necessary, moistened with alcohol) until they are visually clean. Wipe with disinfectant (e.g. disinfectant Bacillol AF from Bode Chemie GmbH (time 30s)) only according to the instructions of use of the respective manufacturer of disinfectant. Only disinfectants with proven effectiveness according to national guidelines should be used. After disinfecting, wipe the examination lamps with a damp cloth to remove potential residue.

Please make sure that the cloth is moistened but not wet, so that no moisture penetrates the openings of the examination lamp.

Make sure that glass and lenses are only cleaned with a dry and clean cloth.

Caution!

The diagnostic instruments are not sterile devices; they cannot be sterilized.

Caution!

The power supply adapter must be unplugged from the socket before cleaning or disinfecting the device!

Never place the examination lamps in liquids! Make sure that no liquids penetrate the housing interior! The article is not approved for machine reprocessing and sterilisation.

This can lead to irreparable damage!

If a reusable device shows signs of material deterioration, it should no longer be reused and should be disposed of/claimed according to the procedures described in the Disposal/Warranty sections

2.2. Accessories

Art. no.: 11340 Transformer 230 V/120 V

Art. no.: 11343 Stand with adapter, including base with 5 castors and tube

Art. no.: 11344 Adapter part between stand and lamp arm

Art. no.: 11345 Wall mount

Art. no.: 11346 Universal clamp

3. Technical specifications

Model: ri-magic®HPLED examination light
Table light/wall light/floor light

Input: 100 V - 240 V/50 - 60 Hz/0.6 A

Output: 5 V/3 A/15 W.

Light source: LED

Classification: type B

Working temperature: 0 °C to +40 °C, relative humidity between 30 and 75%
(non-condensing)

Storage place: -5 °C to +50 °C, up to 85% relative humidity (non-condensing)

Weight: table and wall model: approximately 2.5 kg

Floor model: approximately 7.5 kg

8. Maintenance

The instruments and their accessories require no special maintenance. If an instrument needs to be tested for any reason, please contact Riester, or send it to an authorized Riester dealer in your area, the details of which we will provide you with upon request.

9. Standards

Photobiological safety of lamps and lamp systems DIN EN 62471 IEC 60601-1

10. Disposal of packaging

When disposing of the packaging material, pay attention to the appropriate waste regulations. Keep out of reach of children.

11. Explosion hazard

Do not use this equipment in the presence of flammable anaesthetics, vapours or liquids.

12. Disposal of accessories and device

The service life of this ri-magic® HPLED is 10 years. At the end of its service life, the ri-magic® HPLED, as well as its accessories, must be disposed of in compliance with the guidelines regulating the disposal of such products. If you have questions regarding the disposal of the product, please contact the manufacturer.

13. ELECTROMAGNETIC COMPATIBILITY ACCOMPANYING DOCUMENTS ACCORDING TO IEC 60601-1-2, 2014, Ed. 4.0

During installation and operation of the device, observe the following instructions:

To avoid electromagnetic interference with the operation of the device, do not use the device at the same time as other electronic devices.

To avoid electromagnetic interference with the operation of the device, do not use or stack the device near, on or under other electronic devices.

Do not use the device in the same room as other electronic devices, such as life-supporting devices that have a significant impact on the life of the patient and treatment outcomes, or other measuring instruments or treatment devices that use little electrical power.

Do not use cables or accessories that are not specified for the device, as this can increase the emission of electromagnetic waves from the device and reduce the device's electromagnetic immunity.

Caution:

Medical electrical equipment is subject to special precautions in terms of electromagnetic compatibility (EMC).

Portable and mobile radio frequency communication devices can affect medical electrical equipment. The ME device is intended for operation in a home health care electromagnetic environment and for professional facilities such as industrial areas and hospitals.

The user of the device should ensure that it is operated within such an environment.

Warning:

The ME device may not be stacked, arranged or used directly next to or with other devices. When use close to or stacked with other devices is required, the ME device and the other ME devices must be monitored to ensure intended operation within this configuration. This ME device is intended for use by medical professionals only. This device may cause radio frequency interference or interfere with the operation of nearby devices. It may become necessary to take appropriate corrective measures, such as redirecting or rearranging the ME device or shield.

The ME device assessed does not exhibit any essential performance characteristics in the

sense of EN60601-1, which would present an unacceptable risk to patients, operators or third parties should the power supply fail or malfunction.

Warning: 

Portable RF communications equipment (radios) including accessories, such as antenna cables and external antennas, should not be used in closer proximity than 30 cm (12 inches) to parts and cables of the ri-magic® HPLED specified by the manufacturer. Failure to comply may result in reduced device performance.

Guidelines and manufacturer's declaration – electromagnetic emissions		
The ri-magic HPLED examination lamp is intended for use in an electromagnetic environment as specified below. The customer or user of the headlight should ensure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Groupe 1	The ri-magic® HPLED examination lamp only uses RF energy for its internal functioning. Therefore, its HF emissions are very low and are unlikely to cause interference in nearby electronic devices.
RF emissions CISPR 11	Class A	The ri-magic® HPLED examination lamp is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage network that supplies buildings for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Compliant	

Guidelines and manufacturer's declaration – electromagnetic immunity

The ri-magic® HPLED examination lamp is intended for use in the electromagnetic environment specified below. The customer or user of the ri-magic examination lamp should ensure that it is used in such an environment.

Immunity testing	IEC 60601 test level	Compliance	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	Con: ±8 kV Air: ±2,4,8,15 kV	Con: ±8 kV Air: ±2,4,8,15 kV	Floors should be made of wood, concrete or ceramic tiles. If the floor is covered with synthetic material, the relative humidity should be at least 30%.
Electrically fast transient/burst IEC 61000-4-4	±2 kV for power supply lines	±2kV for power supply lines	The quality of the mains power should be that of a typical business or hospital environment.
Voltage surge IEC 61000-4-5	± 1 kV line to line	±1 kV differential mode	The quality of the mains power should be that of a typical business or hospital environment.
Voltage dips, short interruptions and voltage variations on the input lines for the power supply IEC 61000-4-11	<5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% drop in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5s	<5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% drop in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5s	The quality of the mains power should be that of a typical business or hospital environment. If the user requires needs to continue using the ri-magic HPLED examination lamp during a power interruption, it is recommended to supply the ri-magic HPLED examination lamp with power via an uninterruptible power supply or a battery.
Power frequency magnetic field (50 Hz / 60 Hz) IEC 61000-4-8	30A/m	30A/m	Power frequency magnetic fields should have values characteristic of a typical location in a typical business or hospital environment.

NOTE: UT is the AC mains voltage prior to application of the test level.

Guidelines and manufacturer's declaration – electromagnetic immunity

The ri-magic® HPLED examination lamp is intended for use in the electromagnetic environment specified below. The customer or user of the ri-magic® HPLED examination lamp should ensure that it is used in such an environment.

Immunity testing	IEC 60601 test level	Compliance	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz Outside ISM bands	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the ri-magic® HPLED examination lamp, including cables, than the recommended separation distance calculated using the equation applicable to the frequency of the transmitter. Recommended separation distance 26 MHz to 800 MHz 800 MHz to 2.7 GHz Where P is the maximum output power of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be below the compliance level in each frequency range. Interference may occur in the vicinity of devices marked with the following symbol: 
Emitted RF IEC 61000-4-3	3 V/m 26 MHz to 2.7 GHz	3 Vrms	

NOTE 1: At 26 MHz and 800 MHz, the higher frequency range applies.
 NOTE 2: These guidelines may not apply in all situations. The electromagnetic propagation is affected by absorption and reflection by structures, objects and people.

- a The ISM bands (industrial, scientific and medical) between 150 kHz and 80 MHz are 6.765 MHz to 6.795 MHz, 13.553 MHz to 13.567 MHz, 26.957 MHz to 27.283 MHz, and 40.66 MHz to 40.70 MHz.
- b The compliance levels in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range from 80 MHz to 2.7 GHz are intended to reduce the likelihood that mobile/portable communication devices can cause interference if they are accidentally brought into patient areas. For this reason, an additional factor of 10/3 has been added to the formulas that are used to calculate the recommended separation distance for transmitters in these frequency ranges.
- c Field strengths of fixed transmitters, such as base stations for radio (mobile/cordless phones) and land mobile radios, amateur radio, AM and FM radio broadcasts and television broadcasts, cannot be theoretically predicted with accuracy. To evaluate the electromagnetic environment based on fixed RF transmitters, an electromagnetic site survey should be considered. If the field strength measured at the site where the ri-magic® HPLED examination lamp is used exceeds the RF compliance level specified above, the ri-magic® HPLED examination lamp should be observed to verify normal operation. If abnormal performance is detected, additional measures may be required, such as the reorientation or relocation of the ri-magic® HPLED examination lamp.
 Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.

13.1. Disposal

The used medical device must be disposed of in accordance with current medical practices or local regulations on the disposal of infectious biological medical waste.

Batteries and electrical/electronic devices may not be treated as household waste and must be disposed of in accordance with local regulations.

If you have any questions about the disposal of products, please contact the manufacturer or their representative.

14. Warranty

This product was manufactured to the highest quality standards and subjected to a thorough final inspection before leaving our factory.

We are pleased to issue a warranty of **2 years from the date of purchase** on all defects traceable to material or manufacturing defects. A warranty claim is excluded from cases of improper handling or use.

Any defective parts will be replaced or repaired free of charge within the warranty period.

A warranty claim can only be made if the product is accompanied by this warranty card, which is filled out in full and stamped by the dealer.

Please note that warranty claims must be made within the warranty period.

We are of course happy to charge for checks or repairs after the expiry of the warranty period. We also offer free, no-obligation quotes.

In case of warranty coverage or repair, we ask you to return the RIESTER product with the completed warranty card to the following address.

Rudolf Riester GmbH

Repairs dept. RR

Bruckstr. 31

D-72471 Jungingen, Germany

Serial number or batch number:


Date, stamp and signature of the specialist dealer:

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